

FastPlex™ Triplex 1-Step COVID-19 Detection Kit (RT-PCR, RNA extraction Free)

Instructions for Use

**Catalog # 02.01.1020
(96 Tests/kit)**



PreciGenome LLC

Commercial Entity

Address: 2176 Ringwood Ave. San Jose, CA, 95131, USA

Tel: (001) 408-7084602

Email: info@precigenome.com

Website: www.precigenome.com

Table of Contents

1.	Package Specification.....	3
2.	Intended Use.....	3
3.	Product Overview/Test Principle.....	3
4.	Components Included within the Kits.....	4
5.	Reagent Stability and Transportation.....	4
6.	Components Required But Not Included within the Test.....	5
7.	Warnings and Precautions.....	5
8.	Reagent Storage, Handling, and Stability.....	6
9.	Controls Materials.....	7
10.	Collection, Storage and Shipment of Specimens.....	7
11.	Laboratory Procedure.....	7
12.	Interpretation of Results.....	16
13.	Limitations.....	17
14.	Troubleshooting.....	18
15.	Symbols.....	18
16.	Contact Information and Product Support.....	19

1. Package Specification

96 tests/kit

2. Intended Use

The ***FastPlex™ Triplex 1-Step COVID-19 Detection Kit (RT-PCR, RNA extraction free)*** is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test intended for the qualitative detection of nucleic acid from the SARS-CoV-2 in nasal/oropharyngeal swabs from individuals with signs and symptoms of infection who are suspected of COVID-19. This kit meets the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

Results are for the detection of SARS-CoV-2 RNA. The SARS-CoV-2 RNA is generally detectable in respiratory specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

The ***FastPlex Triplex 1-Step COVID-19 Detection Kit (RT-PCR, RNA extraction free)*** is intended for use by trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and in vitro diagnostic procedures. The ***FastPlex Triplex 1-Step COVID-19 Detection Kit (RT-PCR, RNA extraction free)*** meets the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

3. Product Overview/Test Principle

The ***FastPlex Triplex 1-Step COVID-19 Detection Kit (RT-PCR, RNA extraction free)*** is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. The SARS-CoV-2 primer and probe set(s) is designed to detect RNA from SARS-CoV-2 in respiratory specimens from patients with signs and symptoms of infection who are suspected of COVID-19.

The oligonucleotide primers and probes for specific detection of SARS-CoV-2 are selected from regions of Open Reading Frame 1ab (ORF1ab) and the nucleocapsid gene (N) of the SARS-CoV-2 genome. The kit includes primers/probes that are specific for the ORF1ab gene (probe labeled with FAM) and N gene (probe labeled with HEX) of SARS-CoV-2. In addition, the kit also contains primers and a probe (labeled with CY5) for the human RNase P gene as an endogenous internal control for specimen integrity, nucleic acid isolation, amplification and detection.

RNA isolated and purified from upper and lower respiratory tract specimens is reverse transcribed to cDNA and amplified in a Real-time PCR instrument using one-step Master Mix (SARS-CoV-2 Detection Buffer + SARS-CoV-2 Enzyme Mix). Probes consist of a reporter dye at the 5' end and quenching dye at the 3' end. The fluorescent signals emitted from the reporter dye are absorbed by the quencher. During PCR amplification, probes hybridized to amplified templates are degraded by the Taq DNA polymerase with 5'-3' exonuclease activity, thereby separating the reporter dye and quencher and generating fluorescent signals that increase with each cycle. The PCR instrument automatically draws a real-time amplification curve for each optical channel based on the signal change, and calculates cycle threshold (Ct) values (the point at which fluorescence is detectable above background) that are interpreted by the operator to determine the presence/absence of SARS-CoV-2 RNA.

4. Components Included within the Kits

4.1 Rainamp collection, Transport and processing kit

<i>Item No.</i>	<i>Components</i>	<i>Composition</i>	<i>Quantities</i>	<i>Reactions/Tube</i>
1	Virus Preservation Medium with Releasing Agent	Release Reagent	1.2mL×96	96
2	Swab	Swab	96	96

4.2 SARS-CoV-2 Detection Kit (RT-PCR)

<i>Item No.</i>	<i>Components</i>	<i>Composition</i>	<i>Quantities</i>	<i>Reactions/Tube</i>
1	SARS-CoV-2 Detection Buffer	Contains primers and probes for ORF1ab (FAM nucleocapsid (N gene; HEX) and RNase P (IC; Cy5).	672 µl×1	96
2	SARS-CoV-2 Enzyme Mix	Reverse Transcriptase, DNA Polymerase, RNase Inhibitor	96 µl×1	96
3	SARS-CoV-2 Negative Control	Water	500 µl×1	29
4	SARS-CoV-2 Positive Control	gblocks for ORF1ab, N gene and Human RNA for internal control RNase P gene	500 µl×1	29

5. Reagent Stability and Transportation

The Rainamp collection, Transport and processing kit should be stored at Room Temperature before use.

The FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) (in small box) should be stored at -20 °C in the dark and should be transported in a sealed foam box with ice packs. The kit should be stored at -20°C. Unpacked kits should avoid repeated freeze-thaw cycles.

6. Components Required But Not Included within the Test

Consumables not supplied:

- 1.5 mL DNase-free and RNase-free Eppendorf tube
- 0.2 mL PCR tube or strip
- Various models of pipettes and pipette tips (10 μ L, 200 μ L and 1000 μ L tips with filters)
- Centrifuge (can reach to 12,000 rpm)
- Microcentrifuge
- Desktop vortex mixer
- 0.9% saline
- -20°C cold blocks
- 10% bleach
- DNAZap™ (Ambion, cat. #AM9890)
- Disposable powder-free gloves and surgical gowns

Real-Time PCR Instrument(s):

ABI 7500 Real-Time PCR System

7. Warnings and Precautions

This test meets the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

- For in vitro diagnostic use only (IVD).
- This test is only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens.
- Follow standard precautions. All patient specimens and positive controls should be considered potentially infectious and handled accordingly.
- Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- Handle all specimens as if infectious using safe laboratory procedures.
- Please read the instructions carefully prior to operation.
- Samples must be collected, transported, and stored using the exact procedures and conditions recommended by the swab manufacturer and in this package insert. Improper collection, transport, or storage of specimens may impact the performance of this test.

- False positive and false negative results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology.
- Negative results do not preclude infection with SARS-CoV-2 virus and should not be used as the sole basis for treatment or other patient management decision.
- Separate laboratory areas, dedicated to performing predefined procedures of the assay, are required. a) 1st Area: Preparation Area—Prepare testing reagent: b) 2nd Area: Sample processing—Process the specimen and controls: c) 3rd Area: Amplification Area—PCR conducted.
- All materials used in one area should remain in that area and should not be moved or used in other areas. After the assay procedures, the workbench and lab supplies should be cleaned and disinfected immediately.
- All contents in this package are prepared and validated for the intended testing purpose. Do not replace any of the package contents as this will affect the testing performance of the kit.
- Components contained within a kit are intended to be used together. Do not mix components from different kit lots.
- Filter plugged nuclease free pipette tips are required and should be replaced after the addition of each reagent or sample.
- Centrifuge tubes in the assay should be DNase/RNase-free.

8. Reagent Storage, Handling, and Stability

- Store FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) at -20°C in the dark when not in use.
- Use the reagents in FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) within 30 days once opened.
- Completely thaw the FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) before use; spin briefly before use.
- Do not freeze/thaw cycles for FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) more than 3 times.
- The reagents in FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) should be transported in a sealed foam box with ice packs or add dry ice.

9. Controls Materials

Controls – Positive and Negative Controls provided with the test kit include:

- I) SARS-CoV-2 Positive Control: The SARS-CoV-2 Positive Control consists of a mix of gblocks for ORF1ab, N gene and Human RNA for internal control RNase P gene. The positive control should be positive for the ORF1ab gene (Ct equals or less than 39), the N gene (Ct equals or less than 39) and the RNase P gene targets (Ct equals or less than 39). If the results are not positive, the rRT-PCR run is invalid. Positive and negative are defined based on a cutoff of $Ct \leq 39$. Please refer to Table 2.
- II) SARS-CoV-2 Negative Control is RNase free water. SARS-CoV-2 Negative Control is used to detect any reagent or environmental contaminations. The SARS-CoV-2 Negative Control should be negative for ORF1ab (FAM), N (HEX) and RNase P (Cy5). If SARS-CoV-2 Negative Control shows any positive results, it indicates contamination of reagents or samples. All sample results need to be invalidated and results must not be reported. It is recommended to decontaminate the PCR lab and use a new box of un-opened reagent before repeating sample testing.

10. Collection, Storage and Shipment of Specimens

- Adequate, appropriate specimen collection, storage, and transport are important in order to obtain sensitive and accurate test results. Training in correct specimen collection procedures is highly recommended to assure good quality specimens and results.

- Specimen collection: **Rinsing mouth clean with water before sample collection. To realize direct PCR without RNA extraction, swab specimens should be collected only using virus preservation medium with releasing agent in Rainamp collection, Transport and processing kit.** Place swabs immediately into Virus Preservation Medium with Releasing Agent. If using other sample collection kit, RNA extraction maybe required prior to PCR reaction.

- Specimen Transportation: Specimens must be packaged, shipped, and transported at dry ice, overnight. Specimen Storage: Upon receipt, specimens can be immediately processed or stored at 2-8°C for up to 2 hours after collection. For longer term storage, store specimens at -70°C or lower.

11. Laboratory Procedure

a) Equipment Preparation

Clean and decontaminate all work surfaces, pipettes, centrifuges, ABI 7500 Real Time PCR system and other equipment prior to use. The following decontamination agents may be used: 10% bleach, 70% ethanol, or DNAzap™ or RNase AWAY® to minimize the risk of nucleic acid contamination.

Warning: Do not use bleach when using specimen collection systems containing guanidinium isothiocyanate as a stabilizer as it may react with bleach to release toxic cyanide gas.

b) Preparation of the controls

To avoid contamination, the positive control needs to be prepared in an area separate from the amplification and extraction area.

c) Preparation of rRT-PCR Reactions

- 1) Thaw SARS-CoV-2 enzyme mix on ice. Keep the SARS-CoV-2 enzyme mix on ice or cold block all the time during preparation and use, and store it at -20°C immediately after use.
- 2) Thaw all FastPlex Triplex SARS-CoV-2 Detection Kit (RT-PCR) components at room temperature. Vortex all components and briefly spin to collect all liquid at the bottom of the tube.
- 3) If specimens in Virus Preservation Medium with Releasing Agent were frozen, thaw specimens on ice or a cold block.
- 4) **Vortex specimens in Virus Preservation Medium with Releasing Agent 2 mins** and centrifuge for 1 minute. After centrifugation, place the tubes in the cold rack or on ice.
- 5) Prepare all PCR mix in an area separate from the sample preparation area.
- 6) Determine the number of reactions (N is the number of reactions including samples, positive control, negative control) that will be included in the test.
- 7) In a 1.5 mL microcentrifuge tubes (DNase/RNase free) prepare the PCR mix by adding detection mix and enzyme mix based on Table 4 below. Mix the PCR mix thoroughly by vortex. The remaining reagent must be stored at -20°C immediately.

Table 1. Preparation of PCR mix

Components	Volume [μL]	Final Concentration
SARS-CoV-2 Detection Mix	7μl x (N+1)	1×
SARS-CoV-2 Enzyme Mix	1 μl x (N+1)	1×
Total volume [μL]	8 μl x (N+1)	--

- 8) Centrifuge the PCR mix prepared in step 7 for 1 minute to collect contents at the bottom of the tube, and then place the tube in a cold rack.
- 9) Dispense 8 μL of the PCR mix into a 200 μL centrifuge tube. Be sure not to introduce any foam or bubbles into the tubes when aliquoting PCR reaction Mix. Cover the wells and transfer to the sample processing area.
- 10) Add 17 μL of the **specimens in Virus Preservation Medium with Releasing Agent.** to the wells pre-filled with PCR mix in the following order: SARS-CoV-2 Negative Control, specimen(s), and SARS-CoV-2 Positive Control. Be sure to deposit samples with the pipette directly into the reaction mix in PCR tubes. Cover each well, centrifuge at 2000 rpm for 10 seconds, and place into Applied Biosystems ABI 7500 real-time RT-PCR system and record the exact location of controls and each specimen.


11) Running a PCR amplification on ABI 7500 using 7500 software v2.3:

11.1 . Start ABI 7500 real time PCR system: Turn on the computer connected to the system first, then turn on ABI 7500 real time PCR system.


11.2 Load the instrument: Push the tray door to open it, load the prepared plate containing samples and controls into the plate holder in the instrument. Ensure that the plate is properly aligned in the holder. Close the tray door.

11.3. Set up the experiment run:


11.3.3.1. Double-click ABI 7500 icon (7500 software v2.3) on the desktop. A new window should appear, select Create New Document from the menu.




11.3.3.2. Click **Next Experiment** and a new screen will appear as below.




11.3.3.3. Enter the Experiment Name and Click 7500 (96 Wells), Quantitation-Standard Curve, TaqMan® Reagents, Standard (~2 hours to complete a run) as below.




11.3.3.4. Click Plate Setup, and a new screen will appear as below:




11.3.3.5. In the **Define Targets and Samples Tab**, add Targets and Sample(s) as below.



11.3.3.6. In the **Assign Targets and Samples Tab**, select the well containing the samples and controls in **View Plate Layout**, and then **Select all targets in the Assign target(s) to the selected wells**. And choose **None** in **Select the dye to use as the passive reference**.




11.3.3.7. In the same Tab, select the well in View Plate Layout, and assign the sample name in the Assign sample(s) to the selected wells.




11.3.3.8. Click **Run Method**, set the parameters as follows:

- Stage 1: 50°C for 15 min, 1 cycle;

- Stage 2: 95°C for 1 min, 1 cycle;
- Stage 3: 95°C for 15 sec, 55°C for 45 sec, 45 cycles.
- Sample Volume: 25 µL
- Data Collection at Stage 3, Step 2 (55.0 @ 0:45)



11.3.3.9. Click **Run**, click **Save** to save the document and click **START RUN** to run the evaluation.




d) Data Analysis

See below for step-by-step operation of ABI 7500 using 7500 software v2.3 for Data analysis.

1) After the run is completed, click **Analysis**. Click **Amplification Plot** tab and view and adjust the raw data.

- Click **Analysis Settings**, check the baseline Start Cycle and End Cycle. In the **Start Cycle** window should be “3-15.” The **End cycle** window should be 5-20. Users can adjust the values according to the actual situation. Click **Apply Analysis Settings**.
- In the **Plot Settings** window, **Delta Rn vs Cycle** should be selected in **Plot Type**, **Liner** should be selected in **Graph Type**, and **Target** should be selected in **Color**.
- In the **Option** tab, select one target in **Target** and adjust the threshold just above the curve from NTC (noise) .
- Lastly, be sure to click “Re-analyse” icon to update the analysis.



The screenshot shows the 'Analysis Settings' dialog box for 'Cr Settings' in the 7500 Software v2.3 interface. The dialog is titled 'Analysis Settings for 20200525-RX' and has tabs for 'Cr Settings', 'Elag Settings', and 'Advanced Settings'. The 'Cr Settings' tab is active, showing a table of target settings and options for automatic baseline and threshold.


Target	Threshold	Baseline Start	Baseline End
N	5.459577	AUTO	AUTO
ORF1ab	388.751786	AUTO	AUTO
RNase P	AUTO	AUTO	AUTO

Options for 'Cr Settings for N':
 Use Default Settings
 Automatic Threshold
 Threshold: 5.459577
 Automatic Baseline
 Baseline Start Cycle: 3 End Cycle: 15

Buttons: 'Bvert to Default Analysis Settings', 'Apply Analysis Settings', 'Cancel'.

Summary: Wells: 3 Unknown, 0 Standard, 0 Negative Control. 24 Empty.

2) Click **View Well Table** tab to display the cycle threshold (Ct) values.



12. Interpretation of Results

All test controls should be examined prior to interpretation of results. If the controls are not valid, the results cannot be interpreted. The Ct cutoff value of this kit is set as 39 and the end user is required to review fluorescent curves before final interpretation. All the positive curves should be typical S-shape amplification curves or without plateau for weakly positive samples.

1) Positive and Negative Controls


The positive control and negative control for each run are interpreted as described in Table 2 below.

Table 2. Positive and Negative Control Interpretation.


SARS-CoV-2 Positive Control			SARS-CoV-2 Negative Control			Results	Actions
ORF1ab (FAM)	N (HEX)	IC (CY5)	ORF1ab (FAM)	N (HEX)	IC (CY5)		
+	+	+	-	-	-	Valid	Continue to result interpretation
Any one of them shows negative			Not considered			Invalid	rRT-PCR failed, re-run
Not considered			Any one of them shows positive				Extraction, rRT-PCR contaminated, re-run
Result of (-): Ct value >39 or Undetermined Result of (+): Ct value ≤ 39 If there is contamination for the re-run, please perform decontamination procedures.							

2) Examination and Interpretation of Specimen Results

Assessment of specimen test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If the controls are not valid, the results cannot be interpreted. **Table 3** below describes the results interpretation concerning the use of the controls provided with the test. The Ct cutoff value of this kit is set as 39 and the end user is required to review fluorescent curves before final interpretation. **All positive curves should be typical S-shape amplification curves or without plateau for weakly positive samples as below.**



Line-Type



S-shape

ORF1ab (FAM)	N (HEX)	IC (CY5)	Results
+	+	Not considered	SARS-CoV-2 Positive
+	-		
-	+		
-	-	+	SARS-CoV-2 Negative
-	-	-	Invalid

Result of (-): Ct value >39 or Undetermined
 Result of (+): Ct value ≤ 39
 Invalid Result: There is no typical S-shape amplification curve or Ct > 39 or No Ct detected for ORF1ab gene (FAM), N gene (HEX) and internal control (CY5), indicating that the specimen concentration is too low, or there are interfering substances that inhibit the reaction. If upon retest, the result is invalid again, another fresh sample should be collected and tested.

Table 3. Interpretation of Results based on Controls.

13. Limitations

False positive and false negative results can be caused by poor specimen quality, improper sample collection, improper transportation, improper laboratory processing, or a limitation of the testing technology.

Mutation in the target sequence of SARS-CoV-2 or change in the sequence due to virus evolution may lead to false negative results. Improper reagent storage may lead to false negative results.

Use of this assay is limited to personnel who are trained in the procedure. Failure to follow these instructions may result in erroneous results.

The performance of the *FastPlex Triplex 1-Step COVID-19 detection kit (RT-PCR, RNA extraction free)* was established using oropharyngeal swabs. Nasal swabs, nasopharyngeal, and mid-turbinate nasal swabs. Bronchoalveolar lavage fluid specimens are also considered acceptable specimen types for use with the kit. but performance has not been established.














Unverified interfering substances or PCR inhibitors may lead to false negative or invalid results.

The ORF1ab and N gene primer/probes may detect SARS-coronavirus based on *in silico* analysis.

14. Troubleshooting

Problems	Possible Causes	Action
No fluorescent signal is detected in any samples, including positive control	Error in the preparation of the master mixture	Verify each component and ensure the volumes of reagent dispensed during preparation of the master mixture are correct. Repeat PCR mixture preparation.
	Instrument settings error	Verify the rRT-PCR instrument settings are correct.
If the fluorescent signal is detected in a negative control reaction	Contamination of the extraction/preparation area	Clean surfaces and instruments with aqueous detergents, wash lab coats, and replace test tubes and tips in use.
	PCR tube not properly sealed	Ensure plates are sealed correctly.
If the fluorescent signal does not display the sigmoidal characteristic	Components degraded	Use a new batch.
	Poor quality of specimen carrying interferences	Repeat the test with extracted RNA. Or collect the sample after rinsing with water.
	PCR equipment failure	Repeat the test or contact the equipment supplier

15. Symbols

	The product is used in vitro, please don't swallow it.		European union authorization representative
	Validity		Refer to instruction book
	Warning, please refer to the instruction in the annex		Manufacturer
	Product temperature scope		Catalogue number
	Batch number		Contains sufficient for <n> tests
	Avoid overexposure to sun		Date of manufacture
	The product meets the basic requirements of European in vitro diagnostic medical devices directive 98/79/EC		

16. Contact Information and Product Support



Company Name: PreciGenome LLC

Address: 2176 Ringwood Ave. San Jose, CA, 95131, USA

Tel: (001) 408-7084602

Email: info@precigenome.com

Website: www.precigenome.com