



NanoGenerator®
Flex-S Plus:
High throughput
Formulation
Screening platform
for Nucleic Acid
Encapsulated LNPs



Components of Nucleic Acid Encapsulated LNPs



Lipid Components



Cationic/ionizable lipid



Helper lipid



Cholesterol



PEGylated lipid

Genetic Materials




mRNA



SiRNA




DNA Plasmid





Cas9mRNA + sgRNA

Lipid Components and Functions



Cationic/Ionizable Lipids



- Increase nucleic acid encapsulation rate
- Critical for endosomal escape
- Increase transfection efficiency



Structural Lipids


DOPE –facilitate fusion between LNP membranes and cell membranes.
Higher protein expression level.

DSPC – stabilizing lipid membrane structure, enhance nucleic acid encapsulation efficiency




Cholesterol

- Enhancing membrane fluidity
- Increasing LNP stability



PEGylated Lipids

- 0.5-2.5% molar ratio
- Targeting function
 - Increase LNP stability
 - extend circulation time
 - reducing clearance by blood proteins and macrophages
 - immune responses (anti-PEG antibody)
 - reduce cellular uptake and hinder the escape of nanoparticles from endosomes



Generic Material Optimization



mRNA material

- Synthetic cap analogues and capping enzymes
- Regulatory elements in the 5'-untranslated region (UTR) and the 3'-UTR
- Poly(A) tail stabilizes mRNA and increases protein translation
- Modified nucleosides, decrease innate immune activation and increase translation
- Sequence and/or codon optimization increase translation

mRNA vaccines — a new era in vaccinology.
Pardi, N., Hogan, M., Porter, F. *et al. Nat Rev Drug Discov* 2018 17, 261–279
<https://doi.org/10.1038/nrd.2017.243>

siRNA materials

- 2'-Ribose modification
 - 2'-Ome, 2'-F
 - Increase metabolic stability and reduce degradation
- Phosphorothioate (PS)
 - Terminal backbone stabilization
- RISC loading, 5' phosphate modification
 - Prolonged durability on target silencing
- GNA glycol nucleic acid, reducing off-targeting
- 3' End backbone extra stabilization

RNAi-based drug design: considerations and future directions
Tang, Q., Khvorova, A *Nat Rev Drug Discov* 2024 23, 341–36.
<https://doi.org/10.1038/s41573-024-00912-9>







DNA materials

- Sequences optimization
 - Enhance transgene expression
 - Reduce autoimmunity
 - Strong promoter for expression
- Codon Optimization
 - Increase protein expression level
 - Codon preference
 - Secondary structure of resulted mRNA
 - Avoid restriction enzyme sites
 - GC ~40-60%

DNA-Based Nonviral Gene Therapy—Challenging but Promising
Xiaocai Guan, Yufeng Pei, and Jie Song
Molecular Pharmaceutics 2024 21 (2), 427-453
DOI: 10.1021/acs.molpharmaceut.3c00907

High throughput system for LNP preparation



High throughput system for LNP preparation



	PreciGenome NanoGenerator® Flex-S Plus	Sequential microfluidic single- sample mixing	Robotic Liquid Handler
Mixing Methods	Microfluidic mixing	Microfluidic mixing	Pipette mixing
Synthesis Mode	Multi-sample	Single-sample	Multi-sample
Washing Needed	No	Yes	No
Run Time for 96 samples	1 hour	> 4 hours	40min
Sample volume	100 – 500 µL	400 µL – 2 mL	200 µL
Sample conc. range	Flexible	Flexible	Only low lipid concentration (1-2mM lipid)
LNP size difference compared to scale up production	Similar	Similar	20-25% larger
LNP PDI compared to scale up production	Similar	Similar	20-25% larger
EE% compared to scale up production	Similar	Similar	20-25% less
Protocol optimization	Well developed	Well developed	Intense (ratio, speed, concentration, tip choice, etc.)

Features of Flex-S Plus



- The Flex-S Plus System facilitates the rapid screening of nanoparticle formulations and early-stage mRNA candidates, offering a substantial increase in project efficiency.
- With a max throughput of 32 samples per run, 96 samples per hour, the Flex-S Plus greatly streamlines screening processes. It offers comprehensive automation of complex protocols, enabling users to concentrate on other laboratory duties.
- The system also permits experimentation with as little as 20 μ l of payload reagent (e.g. mRNA) while providing control over collection volumes. This allows users to optimize the use of valuable materials.


NanoGenerator® Flex-S Plus



	NanoGenerator® Flex-S/Flex-S Plus	Syringe Pump Systems	Tubing Connection Systems
Dead volume per sample	< 20 µl	0.5 mL	0.5 - 1 mL
Source of dead volume	Micro-channel in the mixing Chip	Syringe, connector, and/or mixing chip	Tubing, connector, and mixing chip
Typical production volume	100 - 500 µL	1 – 10 mL	1 – 10 mL
Minimum input volume (Aqueous :Lipid = 3:1)	Aqueous : 75ul Lipid: 25ul	Aqueous: 1 mL Lipid: 0.5 mL	Aqueous: 1 mL Lipid: 0.5 mL
Estimated minimum mRNA cost	\$50	\$660	\$660
Estimated minimum mRNA cost	\$50	\$660	\$660




NanoGenerator® Flex-S



NanoGenerator® Flex-S Plus


NanoGenerator® Flex-S Plus



- Rapid screening of LNP formulations
- Rapid screening of mRNA/siRNA
- 32 samples per run
- 96 samples within one hour
- Disposable consumables

Model	Flex-S	Flex-S Plus
Multi-sample per run	1 – 4	(1 – 8) × 4 per run Up to 96 samples per hour
Full automation	N/A	Yes
Library preparation	N/A	Optional
Throughput	0.1 – 0.5 ml per sample	0.1 – 0.5 ml per sample
Total flow rate	3 ml/min, 4 ml/min	3 ml/min, 4 ml/min
Flow rate ratio	3:1, 4:1	3:1, 4:1
Size range	40 – 200 nm	40 – 200 nm
PDI	0.05 – 0.2	0.05 – 0.2
Encapsulation efficiency	Up to 99%	Up to 99%
Payload	DNA, mRNA, siRNA, Protein, small mol ecules, etc.	DNA, mRNA, siRNA, Protein, small mol ecules, etc.
Dimension	320 mm × 400 mm × 210 mm	630 mm × 570 mm × 660 mm
Weight	8.1 kg	50 kg

NanoGenerator® Flex-S Plus for screening



Sample Workflow:



1. Load samples in 96 well plates;
2. Seal the 96 well plate (optional);
3. Put consumables on the deck: Chips, 96 well plates, pipette tips, and Gaskets;
4. Set parameters in the software and run the program;
5. Collect samples in 96 well plate;
6. Discard/Change consumable.

NanoGenerator® Flex-S Plus



- Robust multi-sample synthesis
- Reliable performance
- Consistent results


Model	Flex-S Plus
Aqueous phase	Sodium acetate buffer, 100mM, pH5.2
Solvent phase	Lipidflex, 15mM in ethanol
Parameters	3.3ml/min, FRR 3:1, 200µL



- Lipid formulation screening
- 96 samples < 1hour
- 96-well Plate format

Model	Flex-S Plus
Aqueous phase	RNA in Sodium acetate buffer, 100mM, pH5.2
Solvent phase	Different lipid formulation

NanoGenerator[®] - Nanoparticle Synthesis System



PreciGenome (Confidential)

NanoGenerator® Scaling Up



- Transferable results from early screening (Flex-S/Flex-S Plus, 0.1ml) to pre-clinical development (Pro, 200ml), then commercial production (Max: 1L, MAX 40L/H: >20L)



Flex-S/Flex-S Plus: 0.1 – 0.5 ml per sample




Flex-M: 1 – 12 ml
Flex-M Premium: 1 – 75ml



Pro: 2 – 200 ml






MAX : 50ml – 1L
MAX (40L/H): >20L



Reagents	
Aqueous phase	Sodium acetate buffer (100mM, pH5.2)
Solvent phase	LipidFlex, 15mM in ethanol

NanoGenerator[®] — Scale Up



Reagents	
Aqueous phase	Sodium acetate buffer (100mM, pH5.2)
Payload	RNA (~600 nt)
Solvent phase	LipidFlex RNA-LNP kit

System Benefits

High Throughput & Efficiency



- Multiple sample (1/4/32) per run.
- Runtime <5 min for 4 samples, 48/96 samples per hour.

Regulatory Compliance



- Intuitive software (21 CFR Part 11 compliant)
- Single-use mixing cartridge

Scalable & Reproducible



- Direct transfer from discovery to clinical manufacturing
- Reproducible manufacturing

Automation



- Automated workflow
- Real-time data monitoring & recording
- Electronic batch records

High Yield



- Small reagent volume (minimum 50 μ l) for each sample.
- Save up to 80% of RNA/lipid cost

Custom Design & Service



- On-site 3Q installation & qualification
- Custom design & OEM